

# **Therapy Chair Serie II – a good way to take care**

Actualway Product Portfolio



**FRESENIUS  
MEDICAL CARE**

# Actualway Serie II Therapy Chair

## The manufacturer

Actualway is a Portuguese manufacturer of medical equipment for care & treatment such as Therapy, Sampling/Blood Donation, Gynaecology, Multi-functional Chairs, Couches and many accessories.

Fresenius Medical Care cooperates since 2018 for offering functional, comfortable and economic dialysis treatment chairs with high stakeholder satisfaction in Europe and abroad.

Actualway has been increasing its distribution all over the world thanks to their promotion efforts. Actualway's sales strategy is based on building a long term relationship with all their partners, providing quality and service.

## Actualway certification

Actualway products are manufactured in Portugal, tested according to applicable norms.

Actualway is certified

- QM-certification according to ISO 9001
- QM- certification according to EN ISO 13485

**The Actualway vision:** The best way to take care is to get our partner's and patient's satisfaction meeting their needs.

**The Actualway promise:** We produce reliable, strong and comfortable products, because we do care.

We provide quality, after-sales service, reliable products adapted to the needs of patients and users, and an excellent delivery time. Come and join us, so that we can expand business together!



*Castors with central brake*



*Stainless steel arm rest support*



*Vertical and lateral armrest adjustment*

**Actualway Therapy Chair was developed to give comfort, high level of security and wellness during treatments:**

- Dialysis
- Blood transfusion

**Features and advantages**

- 4 independent motors
- Trendelenburg position (-12°) in 12 sec.
- Adjustable height 55 to 75 cm
- Electrical foot rest adjustment
- Cervical pillow
- Lockable handset
- 4 castors with central brake, 10 cm diameter
- Stainless steel arm rest support
- Push handle
- Equipotential bonding cable

<b>Technical specifications:</b>	
Power supply	100 V–240 V; 50/60 Hz
Protection classification	IPX4
Handset protection	IPX6
Maximum Load capacity	250 kg
<b>Logistic data:</b>	
Net weight	80 kg
Gross weight	110 kg
Pallet	80 x 120 x 160 cm
<b>Dimensions:</b>	
Standard	217 x 54 x 55–75 cm
Total width including arm rests	88 cm
Frame	105 x 65 x 25 cm

**Options**

- Special length and width: 217 x 60 cm
- Memory foam upholsteries
- Reading lamp
- PVC plastic foot protection
- PVC plastic legs protection
- Lumbar pillow
- Rechargeable battery



*Foldable tray*



*IV serum pole (Stainless steel/ Epoxy)*



*Paper roll holder*




*ABS bottom cover*



*Reading lamp*

# Colour options




Snow




Olive




Canard



Myosotis



Sapphire



Sand




Pepper



Brown



Macadam




Mist




Light Grey




Pottery



Cherry



Garnet



Coal

## Upholstery colour

### Actualway Dialysis Chair Serie II

F00008500

#### Actualway – Therapy Chair Serie II

(4 independent motors: height, back rest, leg rest and foot rest adjustment, lockable handset, Trendelenburg position (handset operated), vertical and horizontal armrest adjustment, 4 castors with central locking brake (10 cm diameter), weight capacity: 250 kg + extra accessories, equipotential bonding cable, pillow and foot rest in upholstery colour)

Myosotis

### Actualway Therapy Chair Colour Options (new orders only)

F00008507 **Actualway Serie II ABS bottom cover**

F00008508 **Actualway Serie II IV serum pole**

F00008509 **Actualway Serie II ABS foldable table**

F00008515 **Actualway Serie II reading lamp**

F00008516 **Actualway Serie II paper roll holder**

F00008518 **Actualway Serie II battery**

F00008519 **Actualway Serie II upholstery memory foam**

### Actualway Chair Serie II Instruction for Use (IFU) Please note: If not otherwise indicated dialysis chairs will be delivered with IFU in English language

F00008510 **Actualway Chair Serie II IFU English**

F00008517 **Actualway Chair Serie II IFU Portuguese**

F00008512 **Actualway Chair Serie II IFU French**

F00008513 **Actualway Chair Serie II IFU Spanish**

**Actualway Chair Serie II IFU in other language** (Please indicate language)

# Article numbers

## **Actualway Therapy Chair Colour Options** (new orders only)

F00008527 **Actualway upholstery colour Snow**

F00008528 **Actualway upholstery colour Sand**

F00008529 **Actualway upholstery colour Brown**

F00008530 **Actualway upholstery colour Olive**

F00008531 **Actualway upholstery colour Macadam**

F00008532 **Actualway upholstery colour Pottery**

F00008534 **Actualway upholstery colour Cherry**

F00008538 **Actualway upholstery colour Garnet**

F00008540 **Actualway upholstery colour Myosotis**

F00008541 **Actualway upholstery colour Canard**

F00008544 **Actualway upholstery colour Sapphire**

F00008546 **Actualway upholstery colour Light Grey**

F00008548 **Actualway upholstery colour Pepper**

F00008550 **Actualway upholstery colour Coal**

F00008625 **Actualway upholstery colour Mist**

F00008556 **Actualway pillow colour Snow**

F00008557 **Actualway pillow colour Sand**

F00008558 **Actualway pillow colour Brown**

F00008559 **Actualway pillow colour Olive**

F00008560 **Actualway pillow colour Macadam**

F00008561 **Actualway pillow colour Pottery**

F00008562 **Actualway pillow colour Cherry**

F00008562 **Actualway pillow colour Garnet**

F00008563 **Actualway pillow colour Myosotis**

F00008565 **Actualway pillow colour Canard**

F00008566 **Actualway pillow colour Sapphire**

F00008567 **Actualway pillow colour Light Grey**

F00008568 **Actualway pillow colour Pepper**

F00008569 **Actualway pillow colour Coal**

F00008626 **Actualway pillow colour Mist**

F00008525 **Actualway footrest = pillow colour**

F00008526 **Actualway foam armrest = pillow colour**

**Manufacturer:**

Via Actual – Comércio Internacional, Lda.  
CE Castelo da Maia  
Rua Manuel Assunção Falcão, N°337, Armazém 13  
4475-041 Santa Maria de Avioso – Maia  
Portugal  
Tlf./Phone: +351 229 829 663/4  
Fax: +351 229 829 665  
E-mail: comercial@actualway.pt  
www.actualway.pt



**FRESENIUS  
MEDICAL CARE**

Head office: Fresenius Medical Care Deutschland GmbH · 61346 Bad Homburg v. d. H. · Germany  
Phone: +49 (0) 6172-609-0 · Fax: +49 (0) 6172-609-2191  
www.FreseniusMedicalCare.com